



MENTOR KIT

Plan it Simply



Alpha-Bio Tec presents the **Mentor Kit** - a convenient and beneficial new collection of mentor replicate abutments that simplify and ease dentist restoration procedure, by selecting the most suitable abutment for the specific clinical procedures, saving chair time and enable accurate restorative planning and technician work.


THE ALPHA-BIO TEC MENTOR KIT:

- Designed to ensure ease of use for selection of abutment by the dentists and technicians
- Makes the choice of correct abutment quick, simple and accurate
- Enables optimal selection of abutment to suit each and every unique procedure

PRODUCT CONTENTS:

- Autoclave compatible kit box - light, compact and portable
- Full range of mentor abutments (a total of 23 mentor abutments)
- Spare pockets for extra mentor abutments and tools
- Color and laser markings for easy item identification

MENTOR KIT DISPLAY


























Note: The mentor kit (Ref No.5555) is provided with mentor abutments only. The paraguide system and rotation tool should be ordered separately.

Plan it Simply


The mentor abutments are divided into four color coded groups, for easy identification and selection:

- **Gold** - Straight abutments: 10 mentor abutments.
- **Magenta** - 25° abutments: 4 mentor abutments.
- **Blue** - 15° abutments: 5 mentor abutments.
- **Grey** - Screw-retained and ball-attachment: 4 abutments.


 <p>1 M-TLA (5565) TLA (5030) Straight Titanium Abutment</p>	 <p>7 M-TLSP3 (5568) TLSP3 (5368) Simply Straight Titanium Abutment Cuff H3.0mm</p>	 <p>13 M-EAAS15 (5577) EAAS15 (5410) Short Esthetic Anatomic Abutment 15°</p>	 <p>19 M-TLA25 (5580) TLA25 (5130) Angled Titanium Abutment 25°</p>
 <p>2 M-TLSP (5570)* TLSP (5403), TLSP (5404), TLSP (5405) Slim Titanium Abutment with Short platform</p>	 <p>8 M-TLSP2 (5567) TLSP2 (5367) Simply Straight Titanium Abutment Cuff H2.0mm</p>	 <p>14 M-TLAL15 (5576) TLAL15 (5092) Long Angled Titanium Abutment 15°</p>	 <p>20 M-TCT5.5 (5564) TCT5.5-N (5254) TCT 5.5-N Abutment</p>
 <p>3 M-EAAPSS (5573) EAAPSS (5408) Short Straight Esthetic Anatomic Abutment</p>	 <p>9 M-TLSP1 (5566) TLSP1 (5366) Simply Straight Titanium Abutment Cuff H1.0mm</p>	 <p>15 M-TLA15 (5575) TLA15 (5090) Angled Titanium Abutment 15°</p>	 <p>21 M-TCT3.5 (5563) TCT3.5-N (5252) TCT 3.5-N Abutment</p>
 <p>4 M-EAAPS (5574) EAAPS (5409) Straight Esthetic Anatomic Abutment</p>	 <p>10 M-EOAPSS (5571) EOAPSS (5406) Short Esthetic Omni Abutment with PS</p>	 <p>16 M-EAAH25 (5583) EAAH25 (5415) High Esthetic Anatomic Abutment 25°</p>	 <p>22 M-TCT1.5 (5562) TCT1.5-N (5222) TCT 1.5-N Abutment</p>
 <p>5 M-EOAPS (5572) EOAPS (5407) Esthetic Omni Abutment with PS</p>	 <p>11 M-EAAH15 (5579) EAAH15 (5412) High Esthetic Anatomic Abutment 15°</p>	 <p>17 M-EAA25 (5582) EAA25 (5414) Esthetic Anatomic Abutment 25°</p>	 <p>23 M-TB (5561)** TB0.5 (6260), TB2 (6210), TB3 (6280), TB4 (6220), TB5 (6270), TB6 (6290) Straight Ball Titanium Abutment L0.5mm</p>
 <p>6 M-TLSP4 (5569) TLSP4 (5369) Simply Straight Titanium Abutment Cuff H4.0mm</p>	 <p>12 M-EAA15 (5578) EAA15 (5411) Esthetic Anatomic Abutment 15°</p>	 <p>18 M-EAAS25 (5581) EAAS25 (5413) Short Esthetic Anatomic Abutment 25°</p>	<p>Mentor Abutment (M) Original Abutment</p>

* Contains 3 heights grooves that represents all slim abutments | ** Contains 6 heights grooves that represents all ball abutments


MENTOR KIT WORKING PROTOCOL



Verifying the correct angle and cuff height from a range of possible mentors abutments.




Selecting the most appropriate mentor abutment for the final restorative procedure.



Completing the restoration procedure with the original abutment.



Plan it Simply



ABOUT ALPHA-BIO TEC.

For over 25 years Alpha-Bio Tec has been a leader in developing, manufacturing and marketing implants, prosthetics parts, biomaterials and a variety of dental surgical instrumentation. Alpha-Bio Tec believes in making implantology simple, while manufacturing the highest quality products for the global market and providing customers with the best service possible... Implantology!



www.alpha-bio.net

Alpha-Bio Tec's products are CE-marked in accordance with the Council Directive 93/42/EEC and Amendment 2007/47/EC. Alpha-Bio Tec complies with ISO 13485:2012 and the Canadian Medical Devices Conformity Assessment System (CMDCAS).

Alpha-Bio Tec Ltd.

7 Hatnufa St. P.O.B. 3936
Kiryat Arye, Petach Tikva
49510, Israel

Israel

T+972-3-9291000
F+972-3-9235055
sales@alpha-bio.net

International

T+972-3-9291055
F+972-3-9291010
export@alpha-bio.net

EC REP MEDES LIMITED

5 Beaumont Gate, Shenley Hill,
Radlett, Herts WD7 7AR,
England
Tel / Fax: +44 1923859810